

Determination of % citric acid



Molar masses


- $\text{KHF} = 204.23 \text{ g/mol}$
- Dehydrated citric acid = 192.14 g/mol
- Citric acid monohydrate = 210.14 g/mol

Titration

- What is the purpose of a titration?
- Which apparatus are used for a titration?
- How does a burette work?
- How does a pipette work?
- Which information is known and which information is unknown in a titration?
- What is very important in a titration?

$$n_a / n_b = (C_a V_a / C_b V_b)$$

$$n_{\text{KHF}} / n_{\text{NaOH}} = (CV)_{\text{KHF}} / C_{\text{NaOH}} V_{\text{titrasie syfer}}$$



Concentration of Hydrogen ions compared to distilled water



Examples of solutions at this pH

10,000,000	pH = 0	Battery acid, Strong Hydrofluoric Acid
1,000,000	pH = 1	Hydrochloric acid secreted by stomach lining
100,000	pH = 2	Lemon Juice, Gastric Acid Vineger
10,000	pH = 3	Grapefruit, Orange Juice, Soda
1,000	pH = 4	Tomato Juice Acid rain
100	pH = 5	Soft drinking water Black Coffee
10	pH = 6	Urine Saliva
1	pH = 7	"Pure" water
1/10	pH = 8	Sea water
1/100	pH = 9	Baking soda
1/1,000	pH = 10	Great Salt Lake Milk of Magnesia
1/10,000	pH = 11	Ammonia solution
1/100,000	pH = 12	Soapy water
1/1,000,000	pH = 13	Bleaches Oven cleaner
1/10,000,000	pH = 14	Liquid drain cleaner





Titration apparatus



How to use a pipette pump.





How to test burette taps.



Funnel on top of burette and burette readings



Remove the funnel before taking a reading on the burette.



Move your eyes to the level of the liquid surface and read the value of the bottom of the meniscus with 1/10 of the smallest scale marked on the burette.




Determine the % citric acid in fizzy drinks


- Each student works individually.
- Three titrations are done.
- Time management is very important.

- Standardization of NaOH solution (titration 1)
 - $\text{NaOH(aq)} + \text{KHC}_8\text{H}_4\text{O}_4\text{(aq)} \rightarrow \text{KNaC}_8\text{H}_4\text{O}_4\text{(aq)} + \text{H}_2\text{O(l)}$
- Titration of pure citric acid (titration 2)
 - $\text{C}_6\text{H}_8\text{O}_7 + 3\text{NaOH} \rightarrow \text{C}_6\text{H}_5\text{O}_7\text{Na}_3 + 3\text{H}_2\text{O}$
- Titration of citric acid in fizzy drink (titration 3)
- Indicator for all three titrations is phenolphthalein

Standardization of NaOH solution (titration 1)

- Boil 150 mL fizzy drink and let it stand to cool down.
- NaOH is hygroscopic.
- KHF primary standard
- Weigh 0,15g KHF (weigh twice in two conical flasks)
- Dissolve KHF in approximately 50 mL distilled water
- Titrate against given NaOH
- Repeat the titration
- Record the titration numbers of the TWO titrations.
- Standardize the NaOH







Titration of pure citric acid

(titration 2)

- Wash glassware and rinse with distilled water
- Weigh 0.04g citric acid (twice in TWO flasks)
- Dissolve acid in approximately 50 mL distilled water
- Titrate against standardized NaOH
- Repeat the titration
- Record the titration numbers of TWO titrations.
- Determine the concentration of citric acid.
- Determine the mass of citric acid.
- Determine the percentage citric acid in pure citric acid.





Titration of citric acid in fizzy drinks (titration 3)

- Transfer 25 mL fizzy drink with pipette to conical flask.
- Three drops phenolphthalein in conical flask.
- Titrate fizzy drink against standardized NaOH.
- Repeat the titration.
- Determine average titration number of TWO titrations.
- Determine the concentration of citric acid in fizzy drink.
- Determine the mass of citric acid.
- Determine % citric acid in fizzy drink.

- Wash all apparatus.
- COMPLETE questions and calculations on answersheet.
-

Calculations (1) and (2)

$$n_a / n_b = (C_a V_a / C_b V_b)$$

$$n_{\text{KHF}} / n_{\text{NaOH}} = (CV)_{\text{KHF}} / C_{\text{NaOH}} V_{\text{titration number}}$$

$$n_{\text{KHF}} = (CV)_{\text{KHF}} \quad \text{and} \quad n_{\text{KHF}} = m/M_{\text{KHF}}$$

Calculate C_{NaOH} the unknown

$$n_a / n_b = (C_a V_a / C_b V_b)$$

$$n_{\text{citric}} / n_{\text{NaOH}} = (CV)_{\text{citric}} / C_{\text{NaOH}} V_{\text{titration number}}$$

$$n_{\text{citric}} = (CV)_{\text{citric}} \quad \text{and} \quad n_{\text{citric}} = m/M_{\text{citric}}$$

Calculate $\text{mass}_{\text{citric}}$ the unknown

$$(\text{Mass}_{\text{citric}} / \text{mass}_{\text{weighed}}) \times 100 = \% \text{ purity}$$

Calculations (3)

$$n_a / n_b = (C_a V_a / C_b V_b)$$

$$n_{\text{citric}} / n_{\text{NaOH}} = C_{\text{citric}} V_{\text{citric}} / C_{\text{NaOH}} V_{\text{titration number}}$$

Calculate C_{citric} the unknown

Meaning of the concentration of citric acid:

A mol citric acid in 1000 ml

Therefore X mol citric acid in 100 ml is

$$X = (A \times 100) / 1000 \text{ mol}$$

% acid in fizzy drink

Mass in gram per 100 ml

Change mol in 100 ml to (mass in 100 g) by multiplication with M_{citric}

$$\text{Massa} = X \text{ mol} \times M_{\text{citric}}$$

% citric acid in the fizzy drink = answer received from mass